

U.S.A. Made - Free Customer Support

Bariatric Elevating 6-Way Radiographic Table

medical devices

CORE LABS
RADIOGRAPHIC
I L C A

atlas XL^{RAD}

SMOOTH EXTRA WIDE TABLETOP DESIGN -- AN INDUSTRY MUST



- Hospital Grade Radiographic Table
- Space Saving Compact Design
- Elevating and Four-Way Floating Tabletop
- Durable and Reliable Design
- 765 Lb. Capacities
- Priced for Maximum Investment Return
- Solid Steel Scissor Lift - 3/4" Thick
- Durable Powder Coated Surfaces

- Smooth, Flat Tabletop For Patient Comfort
- Aircraft Aluminum & Steel Construction
- Ball Bearing, Linear Bearing & Bronze Bushing Motion Components
- Custom Hospital Grade Low Voltage Actuator - Expected 60,000 Full Capacity Cycle Life
- Slide Out PC Board For Easy Service / Installation

*Shown with optional patient handles

Can be configured, optionally, to command a CoRE labs motorized OTS to automatically maintain a 40" SID

MADE IN



U. S. A.

O R T H O - H O S P I T A L - C L I N I C

