

# Acuity CDR



RadmediX

RadmediX

CHIROPRACTIC X-RAY REIMAGINED


RadmediX 2510 Lance Drive, Kettering, OH 45409 USA  
Tel: 844.723.6334 | Fax: 513.278.0150 | [info@radmedix.com](mailto:info@radmedix.com)  
© 2020 RadmediX Inc. All rights reserved.

CE FDA

## Acuity CDR, Advanced X-Ray Technology

Radmedix has expanded our proven straight arm line with a new advanced x-ray system designed specifically for chiropractic care, the Acuity CDR. The Acuity CDR is a smaller more compact system which comes with a variety of features such as full software integration, magnetic lock brakes, 40kw generator, LED collimator, and a removeable grid.

This x-ray solution can fit in virtually any x-ray room as well as drastically cut down installation time reducing clinical downtime.



The Acuity CDR software is fully integrated to the x-ray generator and x-ray techniques are seamlessly preselected based on the anatomical region of the exam.

## BUCKY ENCLOSURE



Features powerful magnetic lock brakes from handle control to raise and lower the system effortlessly. Fully counterbalanced, ultra small footprint solution.



## TUBE HEAD



Features a dual laser rotating LED collimator with precise alignment and a fixed 40in. SID system with optional Automatic Exposure Control (AEC).



## HEAD TO TOE POSITIONING



The Acuity CDR can seamlessly raise and lower to accommodate all standing exams from skull to lower extremities. This maintains constant tube and bucky alignment, improving workflow.

